

3D Parkinson's disease model of human iPSC-derived dopaminergic neurons for novel drug candidates screening

C. Chiang¹, A. Nicolas¹, J. Jarazo², K.J. Wilschut¹, R. van Vught¹, H. L. Lanz¹, S. J. Trietsch¹, J. C. Schwamborn², J. Joore & P. Vulto¹




¹MIMETAS BV, J.H. Oortweg 19, NL-2333CH Leiden, The Netherlands
²Luxembourg Centre for Systems Biomedicine, University of Luxembourg, Luxembourg
info@mimetas.com | www.mimetas.com

Overview

Parkinson's disease (PD) is a neurodegenerative disorder which is characterized by motor dysfunction and progressive loss of dopaminergic neurons in the substantia nigra. Important advances made in understanding of pathogenic pathways have not yet been translated into any neuroprotective treatments¹. The high heterogeneity of the disease and the lack of preclinical models that adequately recapitulate PD pathology prohibit successful development of neuroprotective therapies.

Here, we show the development of a 3D PD model based on patient-derived induced pluripotent stem cells. Functional dopaminergic neurons were obtained after 6 weeks of differentiation and can be used to identify novel drug candidates.


OrganoPlate® 3D technology




The OrganoPlate® is a high-throughput platform that combines the most recent advances in 3D cell culture and microfluidics². The OrganoPlate® contains 96 tissue chips and is compatible with standard laboratory equipment. The intermediate human neuroepithelial stem cells (NESCs) were seeded in an extracellular matrix (ECM) in the gel channel (pink) and differentiated towards functional dopaminergic neurons in 6 weeks. The adjacent medium channel (blue) supplies the neurons with nutrients.

3D mature neuronal network


The NESCs, generated from healthy control and patient iPSCs, are maintained in their stem cell phase by activation of the SHH and WNT signaling pathways and express the stem cell markers SOX2 and Nestin³.



NESCs are differentiated towards dopaminergic neurons to model the substantia nigra. Immunostaining confirmed the presence of dopaminergic neurons in a 3D environment (neuron marker, β -III-tubulin, and tyrosine hydroxylase (TH) for dopaminergic neurons).




Neurons express the mature neuron marker Synapsin-I after 6 weeks of differentiation




Neuronal activity by calcium imaging

Neuronal firing causes an increase in intracellular calcium that can be detected as a rise in fluorescent signal using the calcium sensitive Fluo-4 dye.



Spontaneous neuronal activity



Differentiated neurons show spontaneous activity from 4 weeks differentiation.

- Time lapses are recorded at 20Hz on an ImageXpress Micro XLS-C system (Molecular Devices)
- Recordings are processed using ImageJ to extract calcium traces of individual neurons
- Variation in the fluorescent trace represents neuronal activity

Conclusions

- We developed a 3D Parkinson's Disease model for screening of novel candidate neuroprotective agents.
- Differentiated neurons express the dopaminergic marker TH and mature neuronal marker Synapsin-I.
- After four weeks of differentiation the neuronal network shows spontaneous firing.

Literature

- Ref. (1) Schapira AHV et al., Lancet 384: 545–555.
 Ref. (2) S.J.Trietsch et al., Lab on a chip, 2013, 8, 3548-3554.
 Ref. (3) Edinson Lucumi Moreno et al., Lab Chip, 2015, 15, 2419.